PRONTEX

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Integra

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GB **INSTRUCTIONS**

This manual is intended to assist you in the safe and efficient operation of BLOOD PRESSURE MONITOR PRONTEX INTEGRA. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled "TIPS ON TAKING YOUR BLOOD PRESSURE".

Measurement accuracy validated in accordance with the International Protocol 2010 of the European society of Hypertension (ESH). Product recommended by "Dabl Educational Trust Limited".

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
CONTENTS

PART NAMES AND PRODUCT COMPONENTS	20	MAKING A MEASUREMENT	28
GENERAL INFORMATION	21	VIEWING/DELETING OF SAVED READINGS (MEMORY FUNCTION)	30
PRECAUTIONS AND CONTRAINDICATIONS	25	TROUBLESHOOTING	31
TIPS ON TAKING YOUR BLOOD PRESSURE	26	RISK OF ELECTROMAGNETIC INTERFERENCE AND POSSIBLE REMEDIES	33
PREPARATION BEFORE USE	27	WARRANTY	35

PART NAMES AND PRODUCT COMPONENTS

GENERAL INFORMATION

INDICATIONS FOR USE


This product is intended for non-invasive measurement of systolic and diastolic blood pressure, determination of pulse rate and calculation of pulse pressure in adults. Device designed to be used comfortably and easily at home, but it's not suitable for use on infants. This device has two memory channels, so to record separately the measurements of two different users. Please always consult with your doctor or physician before use this product to measure blood pressure in child or persons in pregnancy or under pre-eclamptic condition. The product is not intended for use during professional transport of a patient outside a healthcare facility.

METHOD OF MEASUREMENT

This product employs the oscillometric method for measurement of blood pressure and pulse rate. The cuff is connected to the main unit and wrapped around the arm. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the arm in response to each heart beat. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and displayed on the LCD as a digital value.

MEASUREMENT ON INFLATION


Initially the cuff is rapidly inflated to the prescribed value at the beginning of a measurement, during which the inflation mark is displayed. After the pressure reaches this point, the inflation mark disappears and the cuff is inflated more gradually, during which pulse wave is detected. As pulse wave is detected, heart mark flashes. Measurement ends when value of systolic blood pressure is determined and air is automatically exhausted from the cuff.

PERSONAL MODE

Thanks to this special function, the cuff inflation is adjusted based on the latest three readings stored in the selected memory channel, either [1] or [2], so to shorten the measurement time and to improve comfort.

For example, high diastolic blood pressure values from last three measurements will lead to inflation to higher pressure value than the default value. In the event that this value is too high for your current blood pressure, the cuff will rapidly deflate and measurement will restart. To use the personal mode effectively, chose the memory channel before starting a measurement. In fact, this device remembers the last memory channel that was displayed before its turn off. To choose the memory channel before starting a measurement, display a reading saved in the desired memory channel, turn off the monitor once and then start measurement.

WHO CLASSIFICATION SYMBOL

Measured blood pressure values classified as indicated in the following table (WHO guideline). The higher class must be applied when systolic and diastolic pressures fall into different classes. For specific information on your blood pressure, always contact your physician. Never make any judgment on your own based on the classification/indication supplied by this device.

WHO CLASSIFICATIONS	SYS (mmHg)	DIA (mmHg)
Optimal	<120	<80
Normal	120-129	80-84
High Normal	130-139	85-89
Hypertension (mild)	140-159	90-99
Hypertension (moderate)	160-179	100-109
Hypertension (severe)	>180	>110
Isolated systolic hypertension	≥140	<90


WHO: World Health Organization

PULSE PRESSURE DISPLAY

This product calculates and displays the pulse pressure value at the end of each measuring. While systolic blood pressure keeps increasing with age, diastolic blood pressure tends to start decreasing at around age 50 years. Pulse pressure is value subtracting diastolic value from systolic value and therefore pulse pressure tends to increase with age. Pulse pressure has been considered to be related to stiffness of artery and has been studied as one of circulatory risk factors. It is sometimes said that 45 mmHg of pulse pressure is normal value. Pulse pressure only does not represent arterial sclerosis, however, observation of pulse pressure in long term will be important.

PRESSURE COLUMN INDICATION

The measurement result of systolic and diastolic pressures is indicated with the pressure column at the left side of the display. The bottom of the column shows the diastolic pressure and the top of the column shows the systolic. The ranges of the pressure values for each segment are described in the right diagram. The pressure column is also indicated when viewing the saved readings.


BODY MOTION SYMBOL

Blood pressure value taken while moving or talking, cannot be said to be the correct value because potentially affect from these conditions. This product analyzes pulse wave and displays (👤) when body motion is detected. (👤) indicates the measurement results might be affected by body movement. The symbol appears on LCD at the end of measurement together with the measured values. In this case, it is suggested to consider that the obtained data may not be reliable, and to repeat measurement with the calm condition.

IRREGULAR PULSE RHYTHM SYMBOL

Pulse rhythm can be disturbed from talking, moving or arrhythmias. This product displays (👤) when the difference of shortest and longest time intervals between pulse waves is over 25%, indicating irregular pulse rhythm.

MEMORY FUNCTION

The measured values are automatically saved for later review in either one of two memory channels [M1] or [M2]. These two memory channels can be used to save readings of two persons separately or to save morning readings and evening readings separately. Each memory channel can save up to 60 readings. When the number of stored readings reaches 60 (max capability), the oldest reading will be deleted to record new reading. Saved readings are recalled with measurement date and time when clock function is activated. *The clock needs to be activated before measurement for date and time to be saved with the measured values.

Reading results that caused an Error [Er] are not saved.

CARE AND MAINTENANCE

Keep the product clean. Inspect its cleanliness after use. Check and ensure the monitor is not dirty before use. In case the cleaning is necessary, use only a soft dry cloth. DO NOT USE gasoline, paint thinner, or other strong solvents. Since the cuff may absorb perspiration and other fluids, inspect it for stain and discoloration after each use. When cleaning the cuff, use a synthetic detergent and gently rub the surface. Air dry thoroughly. Make sure fluid never gets in AIR HOSE. Do not machine wash or scrub it. When storing the product, do not place heavy, hot or sharp objects on it. Do not fold the cuff forcibly. To disconnect the cuff from the main unit, do not pull

AIR HOSE but hold AIR PLUG and unplug it. Do not coil AIR TUBE CONNECTOR too tightly. Keep the batteries out of reach of children. Do not try to stretch or twist the cuff, and do not bring it into contact with sharp or pointed objects.

Do not replace the parts, disassemble or modify the monitor as it causes malfunction of the monitor or danger to users. Stop using the monitor immediately and contact your dealer or the manufacturer in case any visible damage is found on the monitor. Remove the batteries if the product is to be stored for an extended period of time.

As precaution under potential problem caused by bad storage, crash, uncorrected use, we suggest for a check of your blood pressure monitor every 2 years. This operation may only be performed by the manufacturer or by firms authorized by the manufacturer.

PRODUCT SPECIFICATIONS

REF:	DS-B33-03
Model:	PRONTEX INTEGRA
Operating Principle:	Oscillometric method
Indicator:	13 digits liquid crystal display
Pressure Indicating Range:	0 to 300 mmHg (cuff pressure)
Measuring Range:	50 to 250 mmHg (systolic), 40 to 180 mmHg (diastolic), 40 to 160 bpm (pulse rate)
Accuracy*:	±3 mmHg (cuff pressure), ±5% of reading (pulse rate)
Inflation:	Automatic inflation
Exhaust:	Automatic exhaust valve
Power Supply:	4 x 1.5V LR6 (size AA) alkaline batteries or AC adaptor ADP-W5 series (not included)
Power Consumption:	3W (max.)
Electric Rating:	DC6V/3W, with four LR6 batteries Adaptor; AC100-240V, 50-60Hz, 0.12A, Unit; DC6V, 500mA, with designated AC adaptor
Memory:	2 channels, store 60 readings in each channel and average last 3 records
Operating Condition:	+10 °C to +40 °C, 15% to 85% RH (noncondensing)
Transportation/Storage Condition:	-20 °C to +60 °C, 15% to 85% RH (noncondensing)
Cuff:	Coverage arm circumference from 22 to 42 cm, applied part BF type
Main Unit:	Weight; Approx. 315 g (without batteries) Size; Approx. 114 x 174 x 70.4 mm (W x D x H)
Protection against harmful ingress of water:	IP20; Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.
Protection against electric shock:	(If used with batteries) internally powered equipment / (If used with designed AC adapter) Class II equipment
Mode of operation:	Type BF applied part Continuous operation
Expected Service Life:	5 years

**Specifications are subject to change without notice due to improvements in performance. The measurement accuracy of the product has been proven according to ISO 81060-2 protocol. In the clinical study, K5 was used for the determination of diastolic pressure values at all auscultatory measurements. The user is cautioned that the system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual. This product is intended for use in the environment with one atmospheric pressure.*

KEY TO SYMBOLS

SYMBOL	DESCRIPTION		
	Mark of conformity with Directive 93/42/EEC and subsequent amendments		
	Consult the user manual		
	Keep dry		
	Storage / Operating temperature		
	Storage / Operating percentage of humidity		
	Device with BF type applied part (cuff)		
	Device with Class of Insulation II (if used with the provided AC/DC power supply - not included with the device)		
START/STOP	ON/OFF		
	Production lot		
	Serial number		
	Product ID code		
SET	Clock/Date setting		
M1 - M2	Memory/Recall/Setting		
IP20	<p>Degree of protection of an electrical device against accidental or intentional penetration of foreign solid bodies, and protection against the penetration of liquids.</p> <table border="1"> <tr> <td>1st DIGIT Penetration of solids Protected against solids having a dimension greater than Ø 12,5mm</td> <td>2nd DIGIT Penetration of liquids NOT PROTECTED</td> </tr> </table>	1st DIGIT Penetration of solids Protected against solids having a dimension greater than Ø 12,5mm	2nd DIGIT Penetration of liquids NOT PROTECTED
1st DIGIT Penetration of solids Protected against solids having a dimension greater than Ø 12,5mm	2nd DIGIT Penetration of liquids NOT PROTECTED		
	Manufacturer: NIHON SEIMITSU SOKKI CO., LTD. - 2508-13 Nakago Shibukawa Gunma, 377-0293 Japan		

This device complies with European Directives/Standards applicable for Safety, EMC and Measurement Accuracy. Specifications are subject to change without notice due to improvements in performance.

DISPOSAL OF SPENT BATTERIES (Directive 2006/66/EC)

This symbol on the product indicates that the batteries should not be treated as household waste. Ensuring that the batteries are correctly disposed of contributes to preventing potential negative consequences for health and the environment, which could otherwise be caused by their inappropriate disposal. Recycling materials helps to preserve our natural resources. Take the spent batteries to collection points for recycling. For more information about collection and recycling of spent batteries or the product, please contact your Municipality, your local waste disposal service or the shop where you purchased the device.

WARNING REGARDING THE CORRECT DISPOSAL OF THE PRODUCT ACCORDING TO THE EUROPEAN DIRECTIVE 2012/19/EU-WEEE:

After use, the product must not be disposed of with household waste. It can be delivered to specific separate collection centres provided by the municipal administration, or returned to the distributor when purchasing a new device of the same type and having the same functions. Disposing of the product separately avoids possible negative effects on the environment and human health and allows the recovery of materials of which it is composed in order to obtain significant savings in energy and resources. The symbol on the data label indicates the separate collection of electrical and electronic equipment. **Attention:** An incorrect disposal of electrical and electronic equipment could lead to fines.

PRECAUTIONS AND CONTRAINDICATIONS

Do not use this product without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding. For specific information on your blood pressure, is recommended to contact always your physician and never make any judgment on your own regarding measurement results.

Do not use this product along with implantable and wearable medical electrical equipment such as pacemaker, defibrillator, or electrocardiographic monitor. This product is not also intended to be used with HF surgical equipment. Do not use this product in an explosive environment such as near flammable anesthetics or inside oxygen chamber. The system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual.

Do not use cuffs or accessories other than those specified by the manufacturer. Otherwise, correct measurement readings cannot be obtained. The medical devi-

ce requires special precautions regarding electromagnetic compatibility and must be installed and used in accordance with the information provided with the accompanying documents. The device must be installed and used away from mobile and portable RF communication devices (mobile phones, microwave ovens, transceivers, etc.) that may interfere with it. Do not apply the cuff over wounded arm, arm under an intravascular access or therapy or an arterio-venous shunt, or arm on the side of a mastectomy. Otherwise injury may be resulted. Make sure that inflation of the cuff is not causing prolonged impairment of blood circulation. Also, be cautious about temporary loss of the functions of any other medical equipment if any monitoring equipment is used on the same limb with the blood pressure measuring cuff. To avoid harmful injury due to interfered blood flow from cuff inflation:

- Make sure that AIR HOSE is not twisted, folded and / or occluded before measurement. Otherwise, cuff inflation may not be conducted properly and results too much prolonged;

- Do not make measurements repeatedly.

To avoid any possibility of accidental strangulation, keep this product away from children and do not drape AIR HOSE or the cable around your neck.

⚠ Warning: Choking from swallowing small parts and batteries by children or pets is possible, please keep small parts and batteries at places where children and pets can't reach. Because the product includes precision parts, avoid extreme variations of temperatures, exposure to moisture, shock, dust, lint, and direct sunlight. Do not drop or strike the product. This product is not water resistant. Do not press the display and buttons with force and avoid placing the monitor with display face down. Do not take out batteries or unplug the AC adaptor when device is turned on/working. Make sure to switch off the device before removing batteries or AC adaptor. Do not touch the output plug of AC adaptor during measurement.

TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure is the force exerted by the heart in pumping the blood through the arteries and the resistance by the veins in this flow.

Blood pressure varies all the time during the day and/or because influenced by mental and physical factors, and is never constant.

In general, blood pressure is highest during the working hours and gradually decreases during the afternoon and evening hours. It is low during sleep and increases at a relatively fast rate after arising from bed. Influence factors and causes of changes in blood pressure values are:

- Body movement
- Mental Tension
- Emotions
- Drinking Alcohol
- Recent Urination or Bowel Movement
- Changes in the surroundings such as movement or noise, etc.
- Conversation
- Nervousness
- Eating
- Smoking
- Room temperature

Blood pressure measured at home tends to be lower than when measured in a hospital, clinic or doctor's office.

This is because people tend to be more anxious in hospitals or health care facilities and more relaxed in their home environment (effect known as "hypertension from

Do not touch a patient while touching the battery, battery terminals of the device or the output plug of the AC adaptor, such as when replacing the batteries or plugging/unplugging the AC adaptor. Do not inflate the cuff when it is not wrapped around your arm.

This product is not designed for self-use by unspecified persons in public areas. Do not apply the cuff on the limb which the intravenous drip infusion is implemented. In case the product is stored in the environment with ambient temperature above 40°C or below 10°C, please leave it for at least 2 hours before taking a measurement.

⚠ Warning: Do not modify this device without the authorisation of the manufacturer. None of the electrical or mechanical parts of the device have been designed to be repaired by the user. Failure to comply with the above can compromise the safety of the device and, therefore, of its use.

white coat"). This is also why it is so important to know the average values of our blood pressure when normally measured at home.

Always let interpret your blood pressure reading from a qualified physician.

Depending on your age, weight and general condition, blood pressure can be slightly different. Consult with your doctor on determining what blood pressure is normal for you.

For a correct measuring:

- Before taking blood pressure, rest for approximately five minutes and take your blood pressure while relaxed in a quiet environment.
- Measure blood pressure using the correct posture and do not move or speak during measurement.
- Avoid exercise, eating, drinking alcohol, smoking and other activities that affect your blood pressure right before a measurement (2 hours before).
- Take your blood pressure at the same time every day.
- The ambient temperature should be approximately 20°C when you take your blood pressure. The lower temperature may cause the higher blood pressure value due to contracted blood vessels.
- Reading may vary slightly, depending on posture, whether sitting or lying.

PREPARATION BEFORE USE


Install batteries and activate the clock before use. Measurement can be made with clock being inactivated. However, activate the clock to save measurement date and time with readings for better review. The clock is inactivated as the batteries are removed and subsequently will need to be set again. The product (and clock) can be also operated with AC adaptor (see relevant section of this manual).

1. INSTALLING BATTERIES

Open the cover of battery compartment, at the bottom of the main unit. Install four Alkaline 1.5V "AA" size batteries into the compartment. Make sure that the polarities correspond to the (+) and (-) marks inside the battery compartment. Batteries should be installed or removed as their (-) ends are pushed against the spring. Otherwise, short circuit or leakage of the batteries may be caused by scratch against the spring or inappropriate contact with the spring.

CAUTION! Do not use rechargeable batteries or the batteries other than the ones with designated type. The specified performance of the product may not be retained, if those batteries are used. Close the battery compartment cover. Do not force the battery cover into position. Insert two tabs to the main unit first. Then close the cover. Close the cover completely.

Battery symbol  will be indicated when battery power is low. It is an indication that the batteries need to be replaced. When replacing batteries, all four batteries need to be replaced with new ones; do not mix new and old batteries or different batteries. When the battery symbol  starts flashing, the batteries are depleted completely and measurement cannot be taken. Please have the batteries replaced. The battery symbol stays flashing for 10 seconds and turns off automatically.

The enclosed batteries are for monitoring, and their life may be shorter than that of commercial batteries. These batteries are not covered by WARRANTY of this device.

2. ADJUST DATE AND TIME AND ACTIVATE THE CLOCK

The unit has a date/time function which you can activate when necessary. This function is deactivated when the unit is delivered.

ACTIVATING THE DATE/TIME FUNCTION


To activate the date/time function, press START/STOP BUTTON and then immediately, while the full display is visible, press MEMORY BUTTON M1 or M2. This function can be deactivated again at any time in the same way.

SETTING THE DATE/TIME

The four-digit number representing the year appears in the display first. Press MEMORY BUTTON M1 to increase the year and MEMORY BUTTON M2 to decrease it. Press the SET button to confirm the selected year and proceed to setting the date. Set the month with MEMORY BUTTONS M1 and M2 and confirm this with the SET button. Set the day, hours and minutes in the same way.

CHANGING THE DATE/TIME

Date/time can be changed (re-adjusted) by pressing and holding SET button for approximately 3 seconds. Regardless of the activation or deactivation of date/time function before the re-adjustment, the date/time function becomes activated after the re-adjustment.


OPERATION WITH AC ADAPTOR

AC CONNECTOR is located at the back of the monitor.

⚠ Use the designated AC adaptor only (ADP-W5 series).

It is suggested to take out the batteries when using the adapter. However, the clock is inactivated as the AC adaptor is unplugged and therefore you may insert the batteries to keep the clock activated.

In such case, please make a few measurement without plugging the adaptor to avoid

battery leak. Do not connect the AC adaptor to any power socket where accessibility is not satisfying or not according to law.

⚠ Protection: this device is double insulated and protected against short circuit and overload by a primary thermal fuse.


Enclosures and Protective Covers: equipment enclosed to protect against contact with live parts, and with parts which can become live (finger, pin, hook test).

⚠ To avoid potentially hazardous conditions, overheating, accidental drops and/or damages, we recommend disconnecting the power supply from the socket immediately after having finished using the device.

Do not leave the power supply connected to the socket for a long time.

MAKING A MEASUREMENT

Insert AIR TUBE CONNECTOR into AIR TUBE CONNECTION marked with blue printing ring before starting with blood pressure measurement.


⚠ Make sure that the air hose is not bent, twisted or pushed against the table.

⚠ The cuff is applicable to arm with circumference between 22 and 42 cm. Before measurement check your arm size is included in this size range.


⚠ Choose your memory channel before measurement to use the personal mode effectively.

1. SIT ON A CHAIR AND APPLY CUFF

Find a chair and a table so where to be comfortably seated with your feet flat on the floor and your back and arm well supported.


The cuff should be applied preferably to the left arm (or in any case always at same arm of previous measurements) and the air hose comes out toward your hand.

⚠ Inaccurate readings may result if cuff is wrapped over heavy clothing or tight shirt that may restrict circulation in your upper arm or if sleeve is rolled up.



Attach the cuff with hook and loop fastener with the lower edge of the cuff approximately 2 to 3 cm above the inner elbow.

The cuff should fit snugly, leaving enough space for two fingers to be put between your arm and the cuff.


! If the cuff is wrapped more tightly or loosely than this, inaccurate blood pressure readings may result.

Adjust the position of the cuff so that AIR HOSE is over the inner part of your arm over the brachial artery.

Press the surface of the cuff to make sure that the hook & loop fastens securely.

2. Let your forearm rest naturally on the table and keep the cuff at the height of your heart during measurement.


! If the cuff is at a lower level (or higher) of the heart, the measurement values tend to increase (or decrease).

Take deep breaths and relax. Do not move or speak while taking your blood pressure.

Do not cross your legs during measurement.


Do not hold your breath.

3. TURN ON DEVICE Press "START/STOP" BUTTON.

Air is exhausted from the cuff. Device begins the automatic cuff inflating after three short beeps sound.


! If you need to stop measuring press START/STOP BUTTON. Device will deflate immediately and quickly the air in the cuff and turns off.

Inflation symbol disappears and inflation becomes slow.

As soon as the unit detects the heartbeat, it displays the corresponding symbol on the display and starts the measurement and identification of blood pressure values. As soon as the detection of blood pressure is ended, device stops the inflation phase and begins to bleed air from the cuff and in the meantime the measurement results are displayed with a sound of a buzzer. At measurement results display, the blood pressure/pulse rate and the pulse pressure are alternately indicated as shown in

the following pictures. For information regarding other symbols that might appear at the measurement, please see section GENERAL INFORMATION of this manual. For information on the indication [Er], in case it appears during the measurement, please see section TROUBLESHOOTING of this manual.

4. STORE THE MEASUREMENT RESULT

Press either one of MEMORY BUTTON M1 or M2. The result will be stored in the selected channel.

VIEWING/DELETING OF SAVED READINGS (MEMORY FUNCTION)

1. PRESS MEMORY BUTTON

Press MEMORY BUTTON M1 to review readings saved in the memory channel 1 and MEMORY BUTTON M2 to review those in the channel 2.

The average of the last three measurements is displayed with the indication of A3.

Note: Average will not be displayed unless there are two or more readings saved. In case only two readings are saved, the average of the last two measurements is also displayed with the indication of A3.

2. PRESS THE SAME MEMORY BUTTON AGAIN. The latest reading saved in the selected channel is displayed. The reading of the blood pressure/pulse rate and the

5. TURN OFF THE DEVICE

Press "START/STOP" BUTTON.

The monitor will be turned off. Even if you do not turn off the monitor, it turns off automatically after 2 minutes. The measurements result is stored in the memory as the device is turned off.


Never repeat measurements in quick succession, in fact the congestion of blood might determine false readings. Let rest the arm for at least 5/10 minutes between measurement and next one.

Do not perform continuous measurements during the day, limit to those indicated by your doctor.

pulse pressure are indicated alternately. The display at the bottom shows memory number, measurement date and time alternately for the blood pressure/pulse rate display, and measurement date and time alternately for the pulse pressure display. The latest reading is displayed as memory "No. 1". The results are recalled from the latest to the oldest with operation of MEMORY BUTTON. The bigger memory number indicates the older number.


Press "START/STOP" BUTTON.

Even if you do not turn off the monitor, it displays the result for approximately 2 minutes and turns off automatically.

SWITCHING MEMORY CHANNELS

Pressing MEMORY BUTTON M2 at display of reading in the channel 1 will switch the display to reading in channel 2 and pressing MEMORY BUTTON M1 again will return to the display of result in the channel 1.

In case a reading, which the measurement was taken with the date/time function deactivated, is recalled when the date/time function is activated, date and time on the display are indicated with --/-- and --:--.

DELETING SAVED READINGS

Deleting individual readings

Press the memory button to display its recorded readings. Then, scroll through the readings by repeatedly pressing the memory button until you reach the one to be deleted. Then hold down memory button for 4 seconds until display begins to flash and continue to hold down button for about another 4 seconds, until reading is deleted and the display is cleared [•••].

Deleting all the readings in a memory channel

As previously described in this manual, select and display the average of readings of the memory channel to be deleted.

Then delete it by holding down the memory button continuously until the display is cleared [•••].

TROUBLESHOOTING

???

Blood pressure is extremely high or low

Measurement is not taken with correct posture or the cuff is wrapped incorrectly.

Reconfirm the measurement procedures.

The reading was effected from moving or talking.

Remain still and quiet during measurement.

Measurement was taken right after exercise or other activities influencing the reading.

Measure again after resting for more than 5 minutes. See TIPS ON TAKING YOUR BLOOD PRESSURE.

???

Measured values vary all the time

Mental and physical condition or measurement condition is influencing blood pressure.

Make measurement under the same condition.

???

Measured value is different from that taken at a hospital

Mental condition such as nervousness at a hospital or relaxed feeling at home is influencing blood pressure.

Take records of blood pressure taken at home and consult with your doctor.

Blood pressure is changing 24 hours a day in relation with mental and physical conditions. Eating, drinking, smoking, exercising or bathing will alter blood pressure.

Blood pressure will be also influenced by mental condition such as feeling tensed or relieved and health condition. It is important to know the tendency of blood pressure by making measurement regularly. Make your own rules including conditions and time to make measurement and take blood pressure every day.

[A] Symbol is displayed again in the course of inflation and deflation and inflation of the cuff are repeated

When the initial pressure was too high for blood pressure to be taken or if you move or strain your arm or hand, the cuff is deflated and the measurement will re-start.

This does not indicate any problem or malfunction of the monitor.

Do not move or speak during measurement.

[Er 300] is displayed

Over pressure: blood pressure could not be taken due to moving or speaking although the cuff was inflated to the maximum pressure. **Do not move or speak during measurement.**

[Er -1] is displayed

Blood pressure could not be taken due to moving or speaking.

Do not move or speak during measurement.

[Er -2] is displayed

The cuff is not securely connected to the monitor. Secure the AIR TUBE CONNECTOR. The cuff is not wrapped appropriately.

Reconfirm how to put the cuff on your arm.

[Weak battery icon] is displayed

Batteries are weak. **Replace all batteries with new ones.**

Nothing is displayed

Batteries are depleted. **Replace all batteries with new ones.**

Batteries are inserted incorrectly. **Reinsert the batteries correctly.**

Battery terminals are not clean. **Clean the terminals with dry cloth.**

AC adaptor is not connected. **Connect AC adaptor.**

Clock is not displayed

The clock is inactivated.

Activate date/time function and adjust date and time.

Measurement date and time are displayed with [·:·:·] and [·:·:·]

Measurement was taken before the clock was activated. **Measurement date and time cannot be saved without clock being activated.**


[0 0] is displayed

“START/STOP” BUTTON was being pressed when batteries are being installed or when AC adaptor is being plugged in.

Turn off the monitor once with “START/STOP” BUTTON.

Cuff end is loosened

Put the cuff with its hook and loop fastener side down. Put the cuff end through the metal ring. Fold back the cuff at the metal ring.

RISK OF ELECTROMAGNETIC INTERFERENCE AND POSSIBLE REMEDIES

This section contains information on device conformity with EN 60601-1-2. The PRONTEX INTEGRA Blood Pressure Monitor is an electromedical device that requires special precautions regarding electromagnetic compatibility and needs to be installed and commissioned according to the information specified in the accompanying documents. Portable and mobile radio communication devices (mobile phones, transceivers, etc.) may interfere with the medical device and should not be used in close

If you cannot get correct measurement with the methods above, contact your dealer. Do not disassemble or tamper with the internal mechanism.

The Manufacturer declares that the information contained in this manual is consistent with the technical and safety specifications of the device to which the manual refers. The stated technical specifications are updated to the date of the publication of this document and only apply to the machine to which they are attached. The Manufacturer reserves the right to make changes or improvements without prior notice in this document.

proximity with, adjacent to or on top of the medical device. If such use is necessary and unavoidable, special precautions should be taken so that the electro-medical device functions properly in its intended operating configuration (for example, constantly and visually checking for the absence of anomalies or malfunctions). The following tables supply information regarding the EMC (Electromagnetic Compatibility) characteristics of the electro-medical device.

Manufacturer's guide and declaration - Electromagnetic Emission

The PRONTEX INTEGRA is for use in the below specified electromagnetic environment. The Customer and/or user of the PRONTEX INTEGRA must ensure that the device is used in such an environment.

Emissions test	Conformity	Guide to the electromagnetic environment
CISPR11 Radiated / Conducted Emissions	Unit 1	The PRONTEX INTEGRA device uses RF energy only for its internal operation. Therefore, its RF emissions are very low and do not cause any interference with any nearby electronic device
CISPR11 Radiated / Conducted Emissions	Class [B]	The PRONTEX INTEGRA device is suitable for use in all environments, including domestic, and those directly connected to the public distribution network providing power to environments used for domestic purposes
Harmonic EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker EN 61000-3-3	Compliant	

Manufacturer's guide and declaration - Electromagnetic Immunity

The PRONTEX INTEGRA Blood Pressure Monitor is for use in the below specified electromagnetic environment. The Customer and/or user of the PRONTEX INTEGRA must ensure that the device is used in such environment.

Immunity test	Test level	Level of conformity	Guide to the electromagnetic environment
Electrostatic discharge (ESD) - EN 61000-4-2	± 6kV in contact ± 8kV in the air	The device does not change its state	The floors should be made of wood, concrete or ceramic. If the floors are covered with synthetic material, the relative humidity should be at least 30%
Fast transients/burst EN 61000-4-4	± 2kV power supply ± 1kV for signal conductors	The device does not change its state	The power supply should be that typical of a commercial or hospital environment
Surge EN 61000-4-5	± 1kV differential mode ± 2kV common mode	The device does not change its state	The power supply should be that typical of a commercial or hospital environment
Voltage dips, short interruptions and voltage variations EN 61000-4-11	<5% U _i , [>95% dip in U _i] for 0.5 cycle 40 % U _i , [60% dip in U _i] for 5 cycles 70 % U _i , [30% dip in U _i] for 25 cycles <5 % U _i , [>95% dip in U _i] for 5 s	--	The power supply should be that typical of a commercial or hospital environment. If the PRONTEX INTEGRA blood pressure monitor user requires that the device operates in continuous mode, we recommend using it with an uninterruptible power supply
Power supply frequency magnetic field EN 61000-4-8	3A/m	The device does not change its state	The power supply should be that typical of a commercial or hospital environment
Note U _i is the power supply voltage value			

Manufacturer's guide and declaration - Electromagnetic Immunity

The PRONTEX INTEGRA Blood Pressure Monitor is for use in the below specified electromagnetic environment. The Customer and/or user of the PRONTEX INTEGRA must ensure that the device is used in such environment.

Immunity test	Level indicated by EN 60601-2	Level of conformity	Electromagnetic environment - guide
Conducted immunity EN 61000-4-6	3Vrms 150kHz to 80MHz (for devices that are not life-supporting)	V _i = 3 V rms	The mobile and portable RF communication devices should not be used closer to any part of the PRONTEX INTEGRA device, including cables, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distances $d = [3.5/V_{i}] \sqrt{P}$ $d = [3.5/E_{i}] \sqrt{P} \text{ from 80 MHz to 800MHz}$ $d = [7/E_{i}] \sqrt{P} \text{ from 800 MHz to 2.5 GHz}$
Radiated Immunity EN 61000-4-3	3V/m 80MHz to 2.5GHz (for devices that are not life-equipment)	E _i = 3 V/m	Where P is the maximum rated output power of the transmitter in Watt (W) according to the transmitter's manufacturer and D is the Recommended separation distance in metres [m]. The field strength from transmitters to fixed RF, as determined by an electromagnetic survey of the site ⁴⁾ , may be less than the level of conformity of each frequency range ⁴⁾ . Interference may occur near the devices marked with the following symbol: 
<p>NOTE 1: The highest frequency range is applied to 80 MHz and 800 MHz.</p> <p>NOTE 2: These guidelines may not apply in all situations. The electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.</p> <p>⁴⁾ The field strengths for fixed transmitters such as the base stations for land RT (mobile phones and cordless) and mobile radio, amateur radio, AM and FM radio and TV broadcasters cannot be theoretically predicted and with accuracy. To establish an electromagnetic environment caused by fixed RF transmitters, an electromagnetic survey of the site should be considered. If the field strength measured where the PRONTEX INTEGRA device is used exceeds the above applicable conformity level, the normal operation of the device itself should be placed under observation. Additional measures, such as a different orientation or position of the device, may be necessary in the case of abnormal performance.</p> <p>⁵⁾ The field strength over a frequency range from 159 kHz to 80 MHz should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile radiocommunication devices and the monitor

The PRONTEX INTEGRA Blood Pressure Monitor is intended to work in an electromagnetic environment in which RF radiated interferences are under control. The PRONTEX INTEGRA device customer or operator can help prevent electromagnetic interference by ensuring a minimum distance between the RF mobile and portable communication devices (transmitters) and the PRONTEX INTEGRA as recommended below, in relation to the maximum output power of the radiocommunication devices.

Maximum rated output power of the transmitter W	Separation distance at transmitter frequency m		
	150 kHz at 80 MHz $d = [3.5/V] \sqrt{P}$	80 MHz at 800 MHz $d = [3.5/E] \sqrt{P}$	800 MHz at 2.5 GHz $d = [7/E] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters with maximum rated output power not stated above, the D recommended separation distance in metres (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum rated output power of the transmitter in Watt (W), according to the transmitter's manufacturer.

NOTE 1: The separation distance for the highest frequency range is applied to 80 MHz and 800 MHz.

NOTE 2: These guidelines may not apply to all situations. The electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

WARRANTY

This device is guaranteed against manufacturing defects for a period of 5 years (except cuff accessory guaranteed for a period of 2 years) from the date of purchasing as shown on your receipt (maintain original copy to have right to warranty services).

During this period, the unit will be repaired or replaced free of charge if the fault is due to faulty design or manufacture.

This warranty does not cover damages or malfunctions caused by use improper or contrary to the instructions of this manual and is not applicable to the commercial batteries supplied at the time of purchasing.

Manufacturer/Distributor cannot be held liable for accidental or indirect damage, if changes have been made to the device, unauthorised repairs and/or technical interventions have been performed or any of its parts have been damaged by accident, improper use and/or abuse.

Any unauthorised intervention, even minimally, on the device voids the warranty immediately and, in any, case does not guarantee correspondence with the technical and safety requirements of the Medical Devices Directive 93/42/EEC (and subsequent changes) and related standards.

Dispositivo medico di classe IIa (93/42/EEC e s.m.i.)
Class IIa medical device (MDD 93/42/EEC and subsequent changes)

NIHON SEIMITSU SOKKI CO., LTD. 2508-13 Nakago Shibukawa Gunma 377-0293 Japan

Nissei Healthcare (UK) Ltd. Henfield, BN5 9SJ UK

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Distributore: Safety S.p.A. Via G. di Vittorio, 17 20813 Bovisio Masciago (MB) - Italy