


3.2.2 Turning on and applying the Pulse Oximeter

Put one of your fingers into the rubber opening of the PULSE O2 Pulse Oximeter with nail facing upward, then release the damp. Press the power button to turn the PULSE O2 Pulse Oximeter on. The oximeter will automatically turn off, if there is no finger in the device for more than 16±2 seconds.



3.2.3 Read data from display screen
Display OLED Screen Model PULSE 02 Description The screen display can scroll through four directions with six different display modes by pressing the power button.



Note:
1. When the battery power is at its lowest level, the battery symbol  will be shown, reminding users to replace the batteries.

2. The plethymogram can be considered accurate if the wave symbol is fluctuating regularly.


7 Clinical SpO2 Accuracy

The below table shows statistic distribution of an invasive controlled desaturation study, which guided by ISO80601-2-61, Annex EE, Guideline for evaluating and documenting SpO2 Accuracy in human subjects. The statistic distribution displayed the accuracy distribution between the range of 70% ~ 100%, which may helpful to user.

Device	Item	100-70%	100 - 90%	90-80%	80 - 70%
PULSE 02	Bs	1.0	0.80	1.26	0.86
	Sres	1.94	0.86	1.61	2.91
	Arms	2.18	1.17	2	2.99

The below is the Bland-Altman graphical plot of samples from invasive controlled desaturation study.

Figure 7.1-Bland & Altman



4 Cleaning and Disinfection

4.1 Cleaning

Switch off the power and take out the batteries before cleaning. Keep the exterior surface of the device clean and free of dust and dirt. Clean the exterior surface (display screen included) of the unit with a soft dry cloth. Use 75% medical alcohol to clean the surface by applying a small amount with a dry cloth to avoid the alcohol entering the device.

4.2 Disinfection

Disinfect the device after each use, if multiple patients use the device. Use 75% medical alcohol to clean the surface that was in contact with the patient.

CAUTION: Do not use strong solvent, e.g., acetone.

CAUTION: Never use an abrasive such as steel wool or metal polish.

CAUTION: Do not allow any liquid to enter the product and do not immerse any part of the device in liquid.

CAUTION: Avoid pouring liquid on the device while cleaning.

CAUTION: Do not leave cleaning solution on the surface of the device.

5 Troubleshooting and Maintenance

5.1 Maintenance

• Replace the batteries in a timely manner, if the battery indication is low.

• Clean the surface of the oximeter before it is used for diagnosis of patients

• Remove the batteries from the battery compartment if the oximeter will not be operated for a long time.

• It is best to store the product where the ambient temperature is between -25°C and 55°C and humidity is 15% to 93%.


• Regular inspection is recommended in order to make sure that no obvious damage is present that may affect the safety and performance of the device.

• Do not expose the device to flammable substances, high or low temperatures or humidity levels outside those in the operation conditions.

5.2 Troubleshooting

Problem	Possible Reason	Resolution
Oxyhemoglobin or heart rate cannot be shown.	1. Finger is not inserted correctly. 2. Patient's perfusion is too low to be measured.	1. Retry by inserting the finger thoroughly. 2. Try a few more times to make sure there is no problem with the product itself. Otherwise seek medical help for exact diagnosis.
Oxyhemoglobin or heart rate is unstable.	1. Finger might not be inserted correctly. 2. Finger is trembling, or patient's body is moving.	1. Retry by inserting the finger thoroughly. 2. Try to help the patient keep calm and still.
Oxyhemoglobin or heart rate is outside of the standard range.	1. Finger might not be inserted correctly. 2. Patient's SPO2 & PR is abnormal.	1. Retry by inserting the finger thoroughly. 2. Seek medical help for further examination.
The Oximeter cannot be turned on.	1. Batteries may need replacement. 2. Batteries might be installed incorrectly. 3. The oximeter might be damaged.	1. Please replace batteries. 2. Please reinstall the batteries. 3. Please contact your point of purchase.
The screen suddenly turns off.	1. The device automatically turns off, if there is no signal detected for more than 16 seconds. 2. Batteries may need replacement.	1. This is a normal process. 2. Replace the batteries correctly.

Figura 7.7.2-Linear regression for SpO2 vs SaO2



8 Disposal

Consider the applicable regulations when disposing of the PULSE O2 Pulse Oximeter and batteries. The pulse oximeter must not be disposed of in the domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the pulse oximeter. **Dispose of the batteries separately from the device.** Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop. Dispose of the batteries separately from the device.

9 Certificate of guarantee

Safety S.p.A. guarantees the PULSE 02 pulse oximeter against any manufacturing defects for a period of 2 years from the date of purchase, provided it is returned to the retailer where it was purchased, in its original packaging and accompanied by the purchase receipt. During this period, the unit will be repaired or replaced free of charge if the fault is attributable to design or manufacturing errors. This warranty does not cover damages or malfunctions resulting from improper use or use contrary to the instructions in this manual. It does not apply to the commercial batteries supplied with the product at the time of purchase. The manufacturer/distributor shall not be held liable for accidental or indirect damages if unauthorized modifications, repairs, or technical interventions have been performed on the device, or if any product components have been incidentally damaged due to improper use and/or abuse. Any unauthorized intervention, even minor, immediately voids this warranty and, in any case, does not ensure compliance with the technical and safety requirements of the Medical Devices Regulation (EU) 2017/745 (MDR) (and subsequent regulations) and related standards.to repair it.

10 Manufacturer's EMC Declaration

1.All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE regarding electromagnetic disturbances for the expected service life.
2.Guidance and manufacturer's declaration - Electromagnetic Emissions and Immunity.

Table 8 - Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions - CISPR 11	Group 1
RF emissions - CISPR 11	Class B
Harmonic emissions	Not application
Voltage fluctuations / flicker emissions - IEC 61000-3-2	Not application

Tabella 2 – Guida e dichiarazione del produttore – Immunità Elettromagnetica		
Test di Immunità	IEC 60601-1-2 - Test level	LCompliance level
Scariche elettrostatiche (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Transitori elettrici/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Not applicable
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV; line(s) to earth: ±2 kV; 100 kHz repetition frequency	Not applicable
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle; At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle; And 70% 25/30 cycles Single phase: at 0, 0% 300 cycles	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150kHz to 80MHz: 3Vrms, 6Vrms (in ISM and amateur radio bands); 80% Am at 1kHz	Not applicable
Radiated RF - IEC61000-4-3	10 V/m; 80 MHz – 2.7 GHz; 80 % AM at 1 kHz	10 V/m; 80 MHz – 2.7 GHz; 80 % AM at 1 kHz
Proximity magnetic fields IEC 61000-4-39	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m	30 kHz: 8A/m 134.2 kHz: 65A/m 13.56 MHz: 7.5A/m

NOTE U_i is the a.c. mains voltage prior to application of the test level.

6 Specification

Device Name	PULSE O2
Dimensionas (L×W×H)	63*36*34 mm
Weight	Approx.50g - 60g (including 2 × AAA battery)
Anti-electric Shock Type	Internally powered equipment
Anti-electric Shock Equipment Degree	Type BF
EMC Type	Group I Class B
Enclosure Degree of ingress protection	IP22
Internal Power:	2×AAA 1.5V alkaline batteries
Power Consumption	Below 45mA
Screen	0.96" OLED
SpO2 Display	35-100%
Pulse rate Display	30-250 BPM
Resolution	SpO2: 1% Pulse rate: 1BPM
Measure Accuracy	SpO2: ±3% (70%-100%); Unspecified (<70%) PR: ±2BPM
Data averaging and other signal processing	8s
Data Update Period	1s
Operating Environment	Temperature: +5°C to +40°C Humidity: 15% to 93% non-condensing Air Pressure: 70Kpa-106Kpa
Storage & Transport Environment	Temperature: -25° to +55° Humidity: 15% to 93% non-condensing Air Pressure: 70Kpa-106Kpa
Service life	24 months

Tabella 3 – Guidance and manufacturer's declaration - Electromagnetic Immunity							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (V)	Distance (m)	Immunity test Level (V/m)
	385	380–390	TETRA 400	Pulse Modulation 18 Hz	1,8	0,3	27
	450	430–470	GMRS 460, FRS 460	FM ±5kHz deviation, 1 kHz sinusoidal	2	0,3	28
	710	704–787	LTE Band 13, 17	Pulse Modulation 217 Hz	0,2	0,3	9
	745						
	780						
	810	800–960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse Modulation 18 Hz	2	0,3	28
	870						
	930						
	1720	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UNITS	Pulse Modulation 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400–2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation 217 Hz	2	0,3	28
	5240	5100–5800	WLAN 802.11 a/n	Pulse Modulation 217 Hz	0,2	0,3	9
	5500						
	5785						

Tabella 4 - Guidance and manufacturer's declaration - electromagnetic Immunity		
TEST FREQUENCY	MODULATION	IMMUNITY TEST LEVEL (A/m)
30 kHz	CW	8
134.2 kHz	Pulse modulation a 2.1 kHz	65 b
13.56 MHz	Pulse modulation a 50 kHz	7.5 b
a) The carrier shall be modulated using a 50% duty cycle square wave signal. b) r.m.s., before modulation is applied.		