
Prime Therm

TERMOMETRO DIGITALE

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
DIGITAL THERMOMETER

Instruction manual

Carefully read this manual before use and keep it with the device.

AVVERTENZE

- △ Read instructions thoroughly before using digital thermometer.
- △ Choking Hazard: Thermometer cap and battery may be fatal if swallowed.
- △ Do not allow children to use this device without parental supervision.
- △ Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- △ Do not place thermometer battery near extreme heat (explosion hazard).
- △ Remove battery from the device when not in operation for a long time.
- △ The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease conditions.
- △ Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- △ Do not bite the thermometer. Doing so may lead to breakage and/or injury.
- △ Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
- △ Equipment must be cleaned and disinfected before and after each use.
- △ Do not force the thermometer into the rectum. Stop insertion and abort the measurement when pain is present. Failure to do so may lead to injury.
- △ Do not use thermometer orally after being used rectally.
- △ For children who are two years old or younger, please do not use the devices orally.
- △ If the unit has been stored at temperatures over 5°C 40°C (<41°F >104°F), leave it in 5°C 40°C (41°F-104°F) ambient temperature for about 15 minutes before using it.
- △ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- △ PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the MANUFACTURER. Otherwise, degradation of this equipment could result.
- △ It is not intended for use in the oxygen rich environment and presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
- △ Do not put the thermometer in direct sunlight or with cotton wool, otherwise the accuracy will be affected.
- △ ME equipment should not be cleaned and disinfected while in use.
- △ Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
- △ Manufacturer/Distributor cannot be held liable for accidental or indirect damage, if changes have been made to the device, unauthorised repairs and/or technical interventions have been performed or any of its parts have been damaged by accident, improper use and/or abuse. Any unauthorised intervention, even minimally, on the device voids the warranty immediately and, in any case, does not guarantee correspondence with the technical and safety requirements of the Reg. UE 2017/745 and subsequent changes and related standards.

INTENDED USE

The digital thermometers are intended to measure the human body temperature in regular mode orally, rectally and under the arm. And the devices are reusable for clinical or home use on people of all ages, including children under 8 years old with adult supervision.

INTENDED USER/INDICATION FOR USE

Patient or healthy people who want to measure their body temperature, and medical workers who want to measure the patient's body temperature.
Patient target group: Patients or healthy people who need or want to have their temperature taken.

PLEASE READ CAREFULLY BEFORE USING

This digital thermometer provides a quick and highly accurate reading of an individual's body temperature. To better understand its functions and to provide years of dependable results, please read all instructions first. This appliance conforms to the following standards:

- ISO 80601-2-56 Medical electrical equipment —Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement,
- EN 60601-1-11 Medical electrical equipment —Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare
- Environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

CONTENTS

PRODUCT ILLUSTRATION


Figure 1

PRECAUTION

*The performance of the device may be degraded should one or more of the following occur:

- Operation outside the manufacturer's stated temperature and humidity range.
 - Storage outside the manufacturer's stated temperature and humidity range.
 - Mechanical shock (for example, drop test) or degraded sensor.
 - Patient temperature is below ambient temperature.
 - Portable and mobile RF communications can affect the device.
- The device needs special pre-cautions regarding EMC according to the EMC information provided in the accompany documents. (EMC), come indicato nella documentazione allegata.
*Do not use the devices in the MR environment.

SYMBOL EXPLANATION

| | | | | | | | |
|--|---|--|--|--|--|--|--|
| | Caution | | MR Unsafe | | Type BF applied part | | Medical Device |
| | Refer to instruction manual booklet | | CE conformity marking 0123 represent Notified Body identification number | | IP22: The first num.2: Protected against solid foreign objects of 12.5 mm and greater. The second num.2: Prevent the shell from dropping water vertically when it is inclined at 15° in all directions, and have no harmful effect on the articles | | |
| | Manufacturer | | European Union Authorized representative | | Importer | | Distributor |
| | Date of manufacture | | Batch Code | | Manufacturer's Catalog Number | | Storage and Transportation Humidity limitation: 15-95%RH |
| | 70kPa ~ 106kPa | | Atmospheric pressure limitation: 70-106 kPa | | Storage and Transportation Temperature Limit: -20~-55°C (-4~-131°F) | | |
| | General symbol for recovery/ recyclable | | Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products. | | | | |

SPECIFICATIONS

| | |
|---------------------------------------|--|
| Type: | Digital Thermomete (not predictive) |
| Measure Range: | 32.0°C – 42.9°C (90.0°F – 109.9°F) |
| Precision: | ±0.1°C(±0.2°F) during 35.5°C-42.0°C(95.9°F-107.6°F) at 18°C-28°C (64.4°F-82.4°F) ambient operating range ±0.2°C(±0.4°F) for other measuring and ambient operating range |
| Operating mode: | Direct Mode |
| Display: | Liquid crystal display, 3 1/2 digits |
| Memory: | For storing the last measured value |
| Battery: | One 1.5 V DC. button battery (size LR41or SR41, UCC 392) |
| Battery life: | Approx. 200hours of continuous operation or 1 year with 3 measurements per day |
| Dimension: | 12.4cm x 1.8cm x 1.1cm (L x W x H) |
| Weight: | Approx. 11 grams including battery |
| Expected service life: | Three years |
| Ambient operating range: | Temperature: 5°C-40°C(41°F-104°F) Relative humidity: 15%-95%RH Atmospheric Pressure : 700hPa ~ 1060hPa |
| Storage and transportation condition: | Temperature: -20°C-55°C(-4°F-131°F) Relative humidity: 15%-95%RH Atmospheric Pressure : 700hPa ~ 1060hPa |
| Ingress Protection Rating: | IP 22 |
| Classification: | Tipo BF |
| Contraindication: | No contraindication |

°C/°F SWITCHABLE

Temperature readings are available in the Celsius or Fahrenheit scale (°C/°F); located in the upper right corner of LCD.) With the unit off, press and hold the On/Off Button for approximately 2 seconds to change the current setting.

DIRECTIONS

1. Press the On/Off Button next to LCD display. A tone will sound as the screen shows, followed by last recorded temperature. After showing the self-test temperature, the thermometer is now in the testing mode.
2. Position thermometer in desired location (mouth, rectum, or armpit).
 - a. **Oral Use:** Place thermometer under tongue as indicated by "∩" position shown in Figure 2. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.
 - b. **Rectal Use:** Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1cm (less than 1/2 ") into rectum.
 - c. **Armpit Use:** Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required.
3. The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. The minimum measurement time until the signaling tone (beep) must be maintained without exception. The measurement continues even after the buzzer notification. **So that in order to achieve better body temperature measurement result, recommend to keep the probe in mouth and rectum about 2 minutes, or in armpit about 5 minutes regardless of the beep sound and at least 30 seconds measurement interval should be maintained.**


Figure 2

*Note: Normally the buzzes are "Bi-Bi-Bi- Bi-"; Alarm beeps more rapidly when temperature reaches 37.8°C (100°F) or higher, and the buzzes are "Bi-Bi-Bi— Bi-Bi-Bi— Bi-Bi-Bi"

4. To prolong battery life, press the On/Off Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 10 minutes.

TROUBLESHOOTING

| ERROR MESSAGE | PROBLEM | SOLUTION |
|---------------|--|--|
| | Temperature taken is lower than 32.0°C(90.0°F) | Turn off, wait one minute and take a new temperature via close contact and sufficient rest. |
| | Temperature taken is higher than 42.9°C(109.9°F) | Turn off, wait one minute and take a new temperature via close contact and sufficient rest. |
| | The system is not functioning properly. | Unload the battery, wait for 1 minute and reposition it. If the message reappears, contact the retailer for service. |
| | Dead battery: Battery icon is flashing, can't be measurable. | Replace the battery. Dispose of batteries separately from device. |

BATTERY REPLACEMENT

1. Replace battery when the symbol in the lower right corner of the LCD display starts flashing.
2. Pull battery cover off as shown in Figure 3.
3. Gently pull out plastic circuit board with battery chamber approximately 1 cm (slightly less than 1/2"). (See Figure 4)
4. Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 1.5V DC button type LR41 or SR41,UCC392, or equivalent. Be sure battery is installed with "+" polarity facing up. (See Figure 5)
5. Slide battery chamber back into place and attach cover.


Figure 3


Figure 4


Figure 5

Dispose of the batteries separately from the device

CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the use instruction, periodic readjustment is not required. However, we recommend checking calibration every two years or whenever clinical accuracy of the thermometer is in question. Set the temperature of the high-precision water bath temperature bath to 33.00°C, 37.00°C, 42.00°C. Turn on the thermometer and insert into the high-precision water bath and then check the laboratory accuracy of thermometer. Repeat the previous measurement process after it is stable for 10 minutes. Please send the complete device to the dealers or manufacturer. The above recommendations do not supersede the legal requirements. The user must always comply with legal requirements for the control of the measurement, functionality, and accuracy of the device which are required by the scope of relevant laws, directives or ordinances where the device is used.

CLEANING AND DISINFECTION

1. Pre-clean the surface of the thermometer probe using potable water or wipe it with a clean, soft cloth dampened with potable water to remove surface dirt;
2. Immerse the thermometer probe completely in high-purity water for at least 1 minute;
3. Gently wipe the thermometer with a dry, clean, soft cloth to remove any residue;
4. Repeat steps 2 and 3 three times until no dirt is visible on visual inspection after cleaning;
5. Disinfection:
 - After use in the armpit, mouth or rectum
 - Method A(High level disinfection): immerse the thermometer probe in 0.55% OPA (Ortho-phthalaldehyde), such as CIDEX OPA, for at least 12 minutes under temperature at 20°C (68°F); Apply before use or after axillary use
 - Method B (intermediate disinfection): Using a clean soft cloth dipped in 70% Isopropanol, wipe the probe 3 times, at least one minute for each time.
6. Repeat steps 2 to 4 to remove disinfectant residue;
7. After cleaning and disinfection is complete, check that the product looks normal, then put the thermometer in its dust cover (if any) and store the thermometer in a cool, dry, clean environment away from high temperatures and corrosive substances.

Precautions:

1. The device must be cleaned and disinfected before and after each use.
2. If oral or rectal measurements must be taken, we strongly recommend advanced disinfection.
3. For the disinfectants or cleaning agents mentioned above, please follow their user manuals.
4. The water temperature used in all the above processes should not exceed 45°C (113°F).
5. Do not use benzene, paint thinner, gasoline or other strong solvents to clean the thermometer.
6. Do not attempt to disinfect the sensing probe (tip) of the thermometer by immersing in 70% Isopropanol, OPA or in hot water (water over 45°C/113°F) for long time.
7. Do not use ultrasound and other automatic or semi-automatic methods to clean the thermometer.
8. Do not use disinfection methods other than the recommended cleaning and disinfection procedures, unapproved other disinfection protocols may cause permanent damage to the equipment.
9. During the service life of the product, it can be reused about 3,000 times according to the above cleaning and disinfection procedures.

LIMITED WARRANTY

This device is guaranteed against manufacturing defects for a period of 2 years from the date of purchasing as shown on your receipt (maintain original copy to have right to warranty services). During this period, the unit will be repaired or replaced free of charge if the fault is due to faulty design or manufacture. This warranty does not cover damages or malfunctions caused by use improper or contrary to the instructions of this manual and is not applicable to the commercial batteries supplied at the time of purchasing.

Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products. Dispose of the battery separately from the device.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard JEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

ELECTROMAGNETIC COMPATIBILITY INFORMATION - ELECTROMAGNETIC EMISSION

| The device is intended for use in the electromagnetic environment specified below. The patient or the user of the device should assure that it is used in such an environment. | |
|--|-----------------|
| Emissions test | Compliance |
| CISPR 11 Conducted emission | Not applicable |
| CISPR 11 Radiated emission | Group 1 Class B |
| Harmonic emissions IEC 61000-3-2 | Not applicable |
| Voltage fluctuations / flicker emissions IEC 61000-3-3 | Not applicable |

ELECTROMAGNETIC COMPATIBILITY INFORMATION - ELECTROMAGNETIC IMMUNITY

| The device is intended for use in the electromagnetic environment specified below. The patient or the user of the device should assure that it is used in such an environment. | | |
|--|---|---|
| Immunity test | IEC 60601 test level | Compliance level |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact 1_2 kV, H kV, 1_3 kV, 1_U5 kV air | ± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air |
| Electrostatic transient J burst IEC 61000-4-4 | ± 2 kV for power supply lines 100 k. Hz repetition frequency ± 1 kV for input/output lines | N/A |
| Surge IEC 61000-4-5 | ± 0.5 kV, ± 1 kV differential mode line-line | N/A |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0% UT (100% dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT (100% dip in UT) for 1 cycle at 0° 70% UT (30% dip in UT) for 25/30 cycle at 0° 0% UT (100% dip in UT) for 250/300 cycle at 0° | N/A |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 N/m, 50/60 Hz | 30 N/m, 50/60Hz |

NOTE: UT is the a. c. mains voltage prior to application of the test level

ELECTROMAGNETIC COMPATIBILITY INFORMATION - ELECTROMAGNETIC IMMUNITY

| The device is intended for use in the electromagnetic environment specified below. The patient or the user of the device should assure that it is used in such an environment. | | |
|--|--|--|
| Immunity test | IEC 60601 test level | Compliance level |
| Conducted RF IEC 61000-4-6 | 3 Vrms da 150 kHz a 80 MHz 6 Vrms da 150 kHz a 80 MHz (outside ISM banda)* | N/A |
| Radiated RF IEC 61000-4-3 | 10 V/m to 80 MHz a 2,7 GHz | 10 V/m |
| Proximity magnetic fields IEC 61000-4-39 | 30 kHz, 8 A/m, CW 134,2 kHz, 65 A/m, Pulse modulation 2,1 kHz 13,56 kHz, 7,5 A/m, Pulse modulation 50 kHz | 30 kHz, 8 A/m, CW 134,2 kHz, 65 A/m, Pulse modulation 2,1 kHz 13,56 kHz, 7,5 A/m, Pulse modulation 50 kHz |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are: 6,765 MHz – 6,795 MHz; 13,553 MHz – 13,567 MHz; 26,957 MHz – 27,283 MHz; 40,66 MHz – 40,70 MHz.

Amateur radio bands between 0,15 MHz to 80 MHz are: 1,8–2,0 MHz; 3,5–4,0 MHz; 5,3–5,4 MHz; 7,0–7,3 MHz; 10,1–10,15 MHz; 14,0–14,2 MHz; 18,07–18,17 MHz; 21,0–21,4 MHz; 24,89–24,99 MHz; 28,0–29,7 MHz; 50,0–54,0 MHz.

b) The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 1/0,3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

d) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN RF WIRELESS COMMUNICATIONS EQUIPMENT

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The patient or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

| Frequency (MHz) | Maximum Power (W) | Distance | EC 60601 Test Level | Compliance Level |
|-----------------|-------------------|----------|---------------------|------------------|
| 385 | 1.8 | 0.3 | 27 | 27 |
| 450 | 2 | 0.3 | 28 | 28 |
| 710 | | | | |
| 745 | 0.2 | 0.3 | 9 | 9 |
| 780 | | | | |
| 810 | | | | |
| 870 | 2 | 0.3 | 28 | 28 |
| 930 | | | | |
| 1720 | | | | |
| 1845 | 2 | 0.3 | 28 | 28 |
| 1970 | | | | |
| 2450 | 2 | 0.3 | 28 | 28 |
| 5240 | | | | |
| 5500 | 0.2 | 0.3 | 9 | 9 |
| 5785 | | | | |

Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARNINGS!

This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.

The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

The operator should maintain basic safety and essential performance in regards to EMC

JOYTECH Healthcare Co., Ltd.
NO.365, Wuzhou Road,
311100 Hangzhou, Zhejiang Province,
PEOPLE'S REPUBLIC OF CHINA
Tel:+86-571-81957767

The product is in compliance with the requirements of MDR(EU)2017/745, "0123" is the identification number of notify body.

Shanghai International Holding Corp. GmbH(Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Safety S.p.A.
Via G. Di Vittorio, 17
20813 Bovisio Masciago (MB) - Italy
www.safety.it - www.prontex.it

SPLIT CHARGE